

# Sound insulation for ultrasonic homogenizers

Sonation UCB series

## Effective - Sophisticated - Flexible



### Advantages

- Reduces the noise of ultrasonic homogenizers by about 75%.
- Suitable for almost all manufacturers and models.
- Large integrated viewing window.
- Height-adjustable and swiveling holder for the sonotrode head.
- Integrated interior lighting.
- Automated ventilation.
- Automatic program end detection with optical signaling.

### Technical data

- External dims: 440x440x910mm (WxDxH)
- Internal dims: 350x370x820mm (WxDxH)
- Weight: Approx. 33kg
- Clamping range of the holder: 43-78mm
- Supply voltage: 12V (via included power supply)